





### images with polarized light

## Visioscope® Palm VP 25

- images with & without reflection of the skin surface for a close look at the skin
  - dry skin
  - pigmentation spots
  - skin ageing/sun damage
  - lines & wrinkles
  - pores
  - impurities (redness, acne, comedones, etc.)
  - hair structure, dandruff and condition of the scalp
  - and many more applications
- with freeze function for the image



Technical Data:

Camera: Dimensions: 202 mm length, largest Ø 41 mm, Illumination: 2 x 6 white LEDs, special polarization filters, Area of view: 5.2 x 4.4 mm (corresponds to a zoom of approx. 13 x), Power supply: via monitor unit, Type: CMOS PAL

Monitor: Dimensions: 152 x 80 x 25 mm, weight including battery and camera: 400 g

Monitor: Dimensions: 152 x 80 x 25 mm, weight including battery and camera: 400 g Power supply: Lithium-Ion rechargeable battery (152P CGA-103450A), 3.6 V/3900 mAh, nominal energy: 14.4 Wh, Rechargeable: with the supplied power supply or by USB at a computer Type: 3.5" QVGA, Resolution: 240 x 320 Pixel, Illumination: LED (approx 1000 h lifetime) Power supply: Input: 90-264 V AC, 47-63 Hz,0.5 A (max.), Output: 5-7 V, 10 W max. Dimensions: 78 x 42 x 30 mm (without adapter), weight: approx.150 g Technical changes may be made without prior notice.



# Flexible Devices for the Skin Analysis

### Moisture, Sebum, Pigmentation

### **Wireless Probes**



- Perfect devices for a targeted recommendation of skin care products
- Quick and easy measurement at the point of sale
- The sebum measurement can also be done on hair and scalp.
- The result is immediately shown on the display.

Technical Data: power supply: battery (2x AA-batteries alkaline); dimensions: length 20.2 cm, Ø 41 mm; weight: 150 g

measurement principles: sebum: photometrical (Scale 0-99) moisture: capacitative (Scale 0-99) pigmentation: photometrical (Scale 0-99)

Technical changes may be made without prior notice.

Easy and flexible determination of skin parameters:

- Hydro Pen HP 10 (moisture)
- Skintone Pen TP 20 (pigmentation)
- Sebu Pen SP 15 (sebum with pads)









combination of the 3 Pens and the Visioscope® Palm VP 25: Pen Station PS 100